


Features

- Fast—clone PCR products in only five minutes
- Versatile—ideal for blunt or stickyend cloning
- Positive selection vector—>99% recombinant colonies
- Economical—no expensive blue/ white screening
- Compatible with all common E. coli strains (TOP10 and XL1-Blue)



PCR cloning

One kit for all DNA fragments

The Thermo Scientific™ CloneJET™ PCR Cloning Kit is an advanced positive selection system for high-efficiency cloning of PCR products generated with any thermostable DNA polymerase. Any other blunt or sticky-end DNA fragment can be cloned as well. The kit is ideal for phosphorylated or non-phosphorylated DNA fragments.

Ligation into the included positive selection vector takes only five minutes, yielding more than 99% recombinant clones.

All common laboratory *E. coli* strains can be directly transformed with the ligation product.









Free sample

Want to try
before you buy?
Contact your sales
representation

Cloning procedure

Blunt-ended PCR products generated with a proofreading enzyme are ligated directly into the cloning vector. 3'-dA tailed PCR products are blunted prior to ligation in 5 minutes using an enzyme provided with the kit.

Selection principle

The kit contains the novel, ready-to-use positive selection cloning vector pJET1.2/blunt. The vector contains a lethal restriction enzyme gene that is disrupted by ligation of a DNA insert into the cloning site. As a result, only bacterial cells with recombinant plasmids are able to form colonies. Recircularized vector molecules lacking an insert express a lethal enzyme, which kills the host *E. coli* cell after transformation.


Vector elements

For convenience in mapping and manipulation of the insert, the pJET1.2/blunt cloning vector contains multiple cloning site flanking the insertion site. In addition, the vector contains a T7 promoter for *in vitro* and *in vivo* transcription of the cloned gene.


Cat #	Description	Size
K1231	CloneJET PCR Cloning Kit	20 rxns
K1232	CloneJET PCR Cloning Kit	40 rxns
Related Products		
K0502/3	Thermo Scientific™ GeneJET™ Plasmid Miniprep Kit	50/250 preps
K0691/2	Thermo Scientific™ GeneJET™ Gel Extraction Kit	50/250 preps
F-534S/L	Thermo Scientific™ Phusion™ Green High Fidelity DNA Polymerase	100/500 U
K2710/1	TransformAid Bacterial Transformation Kit	40 transf

Kit components:

- pJET1.2/blunt Cloning Vector
- T4 DNA Ligase
- 2X Reaction Buffer
- DNA Blunting Enzyme
- pJET1.2 Forward Sequencing Primer
- pJET1.2 Reverse Sequencing Primer
- Control PCR Product
- Water, nuclease-free
- Detailed Protocol



Cloning efficiency of sticky-end PCR products. A 976 bp sticky-end PCR product, generated with Taq DNA Polymerase, was ligated into different cloning vectors according to suppliers, protocols. Ligation mixtures (2 µL) were used to transform E.coli DH10B cells. Transformation efficiency of competent cells was 1×10' transformants per µg supercoiled DNA.



Cloning efficiency of blunt-end PCR products. A 976 bp blunt-end PCR product, generated with Pfu DNA Polymerase, was directly ligated into different positive selection blunt-end PCR cloning vectors according to suppliers, protocols. Ligation mixtures (2 μ L) were used to transform E.coli DH10B cells. Transformation efficiency of competent cells was 1×10^7 transformants per μ g supercoiled DNA.



