


# Material Data Sheet

## S.C.A.T. SafetyCaps 2.0



<b>Product Group</b>	SafetyCaps for HPLC solvent supply, Generation 2.0
<b>Intended use / basic function</b>	Blocking solvent vapours during solvent supply to HPLC / UHPLC instruments, prevention of contamination and underpressure.
<b>Cap / bottle size(s)</b>	GL 45
<b>Max. flow rate (air valve)</b>	≤ 400 ml / min (suitable for preparative applications)
<b>Pore size (air valve)</b>	1 µm
<b>Tubing connector type</b>	One-piece fitting with integrated ferrule
<b>Tubing connector sizes</b>	1.6 / 2.3 / 3.2 / 4.76 / 6.35 mm OD tubing (depending on model)
<b>Number of tubing connectors</b>	1 / 2 / 3 / 4 / 6 (depending on model)



A Body: PTFE premium grade - medical grade		Test method
<b>Color</b>	White	-
<b>Specific gravity</b>	2.130 - 2.190 g/cm <sup>3</sup>	ASTM D 792
<b>Water absorption</b>	0.01 %	ASTM D 570
<b>Flammability</b>	V-0	UL 94
<b>Tensile strength</b>	≥ 25 MPa	ASTM D 4894
<b>Elongation</b>	≥ 280 %	ASTM D 4894
<b>Hardness</b>	≥ 54 Shore D	ASTM D 2240
<b>Ball hardness</b>	≥ 23 MPa	ASTM D 785
<b>Compression strength at 1% deformation</b>	≥ 4 MPa	ASTM D 695
<b>Deformation under load (140 kg/cm<sup>2</sup> for 24 hrs. at 23°C)</b>	10 - 13 %	ASTM D 621
<b>Permanent deformation (after 24 hrs. relaxation at 23°C)</b>	6 - 7.5 %	ASTM D 621
<b>Coefficient of static friction</b>	0.08 - 0.10	ASTM D 1894
<b>Coefficient of dynamic friction</b>	0.06 - 0.08	ASTM D 1894
<b>Wear factor K</b>	2,900	ASTM D 3702
<b>Wear coefficient</b>	20,000 - 25,000 cm <sup>3</sup> min 10-8 kg m h	
<b>Thermal conductivity</b>	0.34 W/m·k	ASTM C 177
<b>Coefficient of linear thermal expansion from 25 to 100 °C</b>	12 - 15 10 <sup>-5</sup> / °C	ASTM D 696
<b>Dielectric strength</b>	≥ 30 kV/mm	ASTM D 149
<b>Volume resistivity</b>	10 <sup>18</sup> Ohm·cm	ASTM D 257
<b>Surface resistivity</b>	10 <sup>17</sup> Ohm	ASTM D 257

**Address**  
S.C.A.T. Europe GmbH  
Waldecker Straße 7  
64546 Mörfelden-Walldorf  
Germany

**Contact**  
Phone: +49-(0) 6105 / 30 55 86 - 0  
Fax: +49-(0) 6105 / 30 55 86 - 99  
E-Mail: info@scat-europe.com  
Web: www.scat-europe.com


**Bank Account**  
Bank: Frankfurter Volksbank  
I.B.A.N.: DE 0550 1900 0060 0151 1678  
B.I.C.: FF VB DE FF  
BLZ: 501 900 00  
No.: 6101 539 146

**Legal Information**  
Registered Office: 64546 Mörfelden-Walldorf  
District Court: Darmstadt  
Comm. Register: HRB 55045  
Managing Director: Peter Rebehn  
VAT-ID-No.: DE 227617033

<b>Statement on suitability for contact with foodstuff</b>	<b>FDA approved US regulation</b>
	Code of federal regulation 21 CFR Ch.1; section 177.1550 Perfluorocarbon Resins of the Food and Drug Administration/US.
	<b>EU Regulation</b>
	EU 1935/2004 - 10/2011 on plastic materials and articles to come in contact with food.

<b>B Screw Cap: Polyphenylensulfide (PPS)</b>		<b>Test method</b>
<b>E-Module (bending test)</b>	13,500 MPa	ISO 178
<b>Bending strength</b>	260 MPa	ISO 178
<b>Strain of the outer fibre (maximum load)</b>	2.2 %	ISO 178
<b>E-Module (tension test)</b>	15,000 MPa	ISO 527
<b>Breaking resistance (tension test)</b>	190 MPa	ISO 527
<b>Breaking elongation (tension test)</b>	1.7 %	ISO 527
<b>Izod impact resistance 1U</b>	50 kJ/m <sup>2</sup>	ISO 180 / 1U
<b>Impact resistance 1eU</b>	55 kJ/m <sup>2</sup>	ISO 179 / 1eU
<b>Impact resistance 1eU -30°C</b>	52 kJ/m <sup>2</sup>	ISO 179 / 1eU
<b>HDT / A (1,8 MPa)</b>	270 °C	ISO 75-1/-2
<b>Melting temperature (DSC)</b>	280 °C	ISO 11357
<b>Density</b>	1.650 kg/m <sup>3</sup>	ISO 1183
<b>Flammability, 0.75 mm</b>	V-0	UL 94
<b>Yellow card available</b>	Yes	-

<b>C Valve Housing: Polypropylene (PP)</b>		<b>Test method</b>
<b>Melt flow rate (MFR) at 230 °C and 2.16 kg</b>	10.5 g/10 min	ASTM D 1238
<b>Density</b>	905 kg/m <sup>3</sup>	ASTM D 792
<b>Molecular weight distribution</b>	Broad	-
<b>Isotacticity</b>	Medium	-
<b>Anti block agent</b>	No	-
<b>Slip agent</b>	No	-
<b>Tensile test: stress at yield</b>	36 MPa	ASTM D 638
<b>Tensile test: strain at yield</b>	9 %	ASTM D 638
<b>Tensile test: secant modulus at 1% elongation</b>	1,800 MPa	ASTM D 638
<b>Vicat softening temperature at 10 N (VST/A)</b>	154 °C	ASTM D 1525
<b>Vicat softening temperature at 50 N (VST/B)</b>	95 °C	ASTM D 1525




**Address**  
S.C.A.T. Europe GmbH  
Waldecker Straße 7  
64546 Mörfelden-Walldorf  
Germany

**Contact**  
Phone: +49-(0) 6105 / 30 55 86 - 0  
Fax: +49-(0) 6105 / 30 55 86 - 99  
E-Mail: info@scat-europe.com  
Web: www.scat-europe.com

**Bank Account**  
Bank: Frankfurter Volksbank  
I.B.A.N.: DE 0550 1900 0060 0151 1678  
B.I.C.: FF VB DE FF  
BLZ: 501 900 00  
No.: 6101 539 146

**Legal Information**  
Registered Office: 64546 Mörfelden-Walldorf  
District Court: Darmstadt  
Comm. Register: HRB 55045  
Managing Director: Peter Rebehn  
VAT-ID-No.: DE 227617033

D Valve Membrane: PTFE, virginal, porous	
Material characteristics	<ul style="list-style-type: none"> <li>• Symmetric, porous structure</li> <li>• Hydrophobic on both sides</li> <li>• Free of PFOA</li> <li>• Free of heavy metals</li> <li>• Free of animal substances</li> </ul>
Colour	white
Pore size	1 µm
Typical water inlet pressure (WIP)	> 0.4 bar
Temperature resistance	-200 to 250 °C
Chemical resistance	Almost universal
Conformity	<ul style="list-style-type: none"> <li>• Protection class IP68</li> <li>• ISO 16750-5:2010</li> <li>• Directive 2000/53/EG</li> <li>• Directive 10/2011/EG on plastic materials and articles to come in contact with food (Implementation in Germany by BedGgstV), supplemented by directives 2011/8/EG, 2007/19/EG and 2002/72/EG</li> <li>• Directive 2002/95/EG on the restriction of the use of certain hazardous substances in electrical and electronic equipment</li> <li>• Directive 2003/11/EG relating to restrictions on the marketing and use of certain dangerous substances and preparations</li> </ul>



E Valve Flap: Silicone VMQ	
Colour	white
Hardness (Shore A)	50
Max. flow rate	400 ml / min
Sterilisation resistance	ETO, Gamma, Steam
Cv	0.08
Max. back pressure	300 kPa
Opening pressure	0.5 kPa
Conformity	ISO 109993, USP class VI, BGW

F Tube Fitting: Teflon™ PFA 340		Test method
Melt flow rate	14 g/10 min	ASTM D 3307 / ISO 12086
Melting point	305 °C	ASTM D 4591
Specific gravity	2.15	ASTM D 792
Critical shear rate, 372 °C	50 1/s	-
Tensile strength: 23 °C	25 MPa	ASTM D 3307 / ISO 12086
Tensile strength: 250 °C	12 MPa	ASTM D 3307 / ISO 12086
Ultimate elongation: 23 °C	300 %	ASTM D 3307 / ISO 12086

**Address**  
 S.C.A.T. Europe GmbH  
 Waldecker Straße 7  
 64546 Mörfelden-Walldorf  
 Germany


**Contact**  
 Phone: +49-(0) 6105 / 30 55 86 - 0  
 Fax: +49-(0) 6105 / 30 55 86 - 99  
 E-Mail: info@scat-europe.com  
 Web: www.scat-europe.com

**Bank Account**  
 Bank: Frankfurter Volksbank  
 I.B.A.N.: DE 0550 1900 0060 0151 1678  
 B.I.C.: FF VB DE FF  
 BLZ: 501 900 00  
 No.: 6101 539 146

**Legal Information**  
 Registered Office: 64546 Mörfelden-Walldorf  
 District Court: Darmstadt  
 Comm. Register: HRB 55045  
 Managing Director: Peter Rebehn  
 VAT-ID-No.: DE 227617033

Ultimate elongation: 250 °C	480 %	ASTM D 3307 / ISO 12086
Flexural modulus: 23 °C	590 MPa	ASTM D 790 / ISO 178
Flexural modulus: 250 °C	55 MPa	ASTM D 790 / ISO 178
MIT folding endurance (0,20 mm, 8 mil film)	15,000 cycles	ASTM D 2176
Hardness durometer	D55	ASTM D 2240 / ISO 868
Dielectric strength, short time, 0,25 mm (0,010 in)	80 kV/mm	ASTM D 149 / IEC 243
Dielectric constant, 1 MHz (10 <sup>6</sup> Hz)	2,03	ASTM D 150 / IEC 250
Dissipation factor, 1 MHz (10 <sup>6</sup> Hz)	< 0.0002	ASTM D 150 / IEC 250
Volume resistivity	10 <sup>18</sup> Ohm·cm	ASTM D 257 / ISO 1325
Water absorption, 24 hr	< 0,03 %	ASTM D 570
Weather and chemical resistance	Outstanding	-
Limiting oxygen index	> 95	ASTM D 2863 / ISO 4589
Continuous service temperature	260 °C	-
Flammability	V-0	UL 94

G Cap Seal: VITON (FPM)		Test method
Colour	Black	-
Hardness	73 - 77 Shore A	ASTM D 2240
Operating temperature (air)	-20 to 200 °C	-
Density	1.85 - 1.98 g/cm <sup>3</sup>	-
Tensile strength	12 - 14.86 MPa	ASTM D 412
Elongation at break	175 - 216 %	ASTM D 412
Compression strain 22 h / 200 °C	11 - 12 %	ASTM D 395



**Address**  
 S.C.A.T. Europe GmbH  
 Waldecker Straße 7  
 64546 Mörfelden-Walldorf  
 Germany

**Contact**  
 Phone: +49-(0) 6105 / 30 55 86 - 0  
 Fax: +49-(0) 6105 / 30 55 86 - 99  
 E-Mail: info@scat-europe.com  
 Web: www.scat-europe.com

**Bank Account**  
 Bank: Frankfurter Volksbank  
 I.B.A.N.: DE 0550 1900 0060 0151 1678  
 B.I.C.: FF VB DE FF  
 BLZ: 501 900 00  
 No.: 6101 539 146

**Legal Information**  
 Registered Office: 64546 Mörfelden-Walldorf  
 District Court: Darmstadt  
 Comm. Register: HRB 55045  
 Managing Director: Peter Rebehn  
 VAT-ID-No.: DE 227617033