



Clean print.  
Lean process.

**Epredia PrintMate AS Tissue Cassette Printer**



### Paired for Performance

Eprelia PrintMate tissue cassettes feature an all-new surface finish and reinforced print face design engineered to optimize print quality when used with PrintMate AS cassette printers.

They're offered with PrintMate AS-compatible tubes that are precisely engineered to reduce jamming and prevent downtime



Maximize resources and enhance your chain of custody – streamline your labeling process with the Eprelia™ PrintMate™ AS cassette printer.

### Lean your workflow and avoid errors within the laboratory

- Save time – prints each cassette in ten seconds or less
- Compact footprint allows the printer to fit easily next to the grossing station
- Print only the cassettes you need, right at the point of use
- Eliminate errors endemic to handwriting processes
- A variety of configurations are available to fit your laboratory's needs

### The information you need is right where you want it

- Enjoy the benefits of crisp, clear identifiers printed indelibly on your tissue cassettes
- Maximize your efficiency – printing 2D data matrix barcodes on cassettes allows for the automation of sample identification in downstream processes
- Barcodes printed on cassettes provide a seamless link to LIS records with complete case identification information
- Barcodes can be scanned at the microtome to automate slide printing



The Eprelia PrintMate AS cassette printer can be used for batch printing cassettes with optional automated collection system as shown.

## EpreDia PrintMateAS Specifications (Based on PrintMate AS 900)

Dimensions (D x W x H) without cassette tubes	15 x 17 x 11.4 in 38 x 43.5 x 29 cm
Dimensions (D x W x H) with cassette tubes	15 x 17 x 32 in 38 x 43.5 x 81 cm
Weight without cassette tubes	31 lb, 14 kg
Weight with cassette tubes	37.5 lb, 17 kg
Power requirements	100-240 VAC, 50/60 Hz, 150 VA (512 BTU/h)
Print media type:	Thermal transfer ribbon
Print media capacity:	Approximately 20,000 cassettes per roll
Print speed:	Approximately 10 seconds per cassette
Connectivity:	USB
Print media color:	Black
Cassette Capacity:	
PrintMate AS 150	150 cassettes
PrintMate AS 450	450 cassettes
PrintMate AS 900	900 cassettes

## Ordering Information

	QTY	Order Number
<b>EpreDia Printmate AS Cassette Printer</b>		
Includes PrintMate AS, USB cable, power cable, 2, 6 or 12 cassette hoppers, one roll of tape (installed), EpreDia™ LabWriter software™ CD, operator guide		
PrintMate AS 150 Cassette Printer		A83800001
PrintMate AS 450 Cassette Printer		A83800002
PrintMate AS 900 Cassette Printer		A83800003
<b>EpreDia Printmate AS Cassette Printing System</b>		
Includes PrintMate AS, USB cable, power cable, 2, 6 or 12 cassette hoppers, one roll of tape (installed), EpreDia LabWriter software CD, operator guide, all-in-one PC, digital imager		
PrintMate AS 150 Cassette Printing System		A83800004
PrintMate AS 450 Cassette Printing System		A83800005
PrintMate AS 900 Cassette Printing System		A83800006

## Accessories

Collection system (includes cable, 4 trays for cassettes, output guides)	A79510014
Collection trays (box of 4)	A79510022
Thermal transfer ribbon, box of 5 rolls	A79510013
All-in-one touch screen PC	B81300033
Skanner Kit & Packaging	B35310050
Cassette hopper (ea)	A79510024

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk  
Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it  
Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl  
Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es  
Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2023 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at fishersci.com/trademarks.

